

PRL PH PAPER ROLL



**\$33.90 inc GST**

pH Paper Roll - pH 5.5 - 8.0

15 Feet per roll

**SKU:** 500-354 | **Categories:** [Brand](#), [PRL](#) | **Tags:** [paper](#), [pH](#), [strips](#), [testing](#)

PRODUCT DESCRIPTION

pH Paper Roll - pH 5.5 - 8.0

**ADDITIONAL INFORMATION**

<b>Weight</b>	40 g
<b>Dimensions</b>	55 × 65 × 15 mm
<b>Size</b>	15 Ft
<b>Indications</b>	Practitioner tools

**Usage:**

- Dip into solution
- Remove at once
- Immediately compare to color chart

Nutritional supplements should be stored in a cool, dry place, with caps tightly closed, and out of reach from children. Take only as directed.

**THERE ARE NO REVIEWS YET.**

